

# ReVel™ ventilator accessory catalog



CareFusion

# Table of contents

ReVel ventilator .....	2
PalmTop Monitor (PTM™) II graphics monitor .....	2
PTV docking stations, stands and accessories...	3–4
PTV power accessories and batteries .....	4
PTV miscellaneous parts and accessories.....	5–6
PTV ReVel circuits .....	7–8



## ReVel™ ventilator accessories



The ReVel™ ventilator, part of the PalmTop™ Ventilator (PTV®) family, provides portable ventilation on the fly with critical tools and performance for your pediatric ( $\geq 5$  kg) to adult patients. The patented ActiveCore™ gas delivery system with unprecedented miniaturization allows a single ventilator to deliver critical care ventilation at the golden hour, on the scene, en route and through all points of care.

### ReVel ventilator



**19260-001**

#### ReVel ventilator

ReVel ventilator, main battery, power adapter w/power cord, oxygen (O<sub>2</sub>) hose, carrying case, O<sub>2</sub> sensor cable, O<sub>2</sub> sensor, O<sub>2</sub> tee and operator's manual

### PalmTop Monitor (PTM™) II graphics monitor



**19400-001**

#### PTM graphics monitor



**99260-001**

#### ReVel system w/floor stand and docking station

ReVel ventilator, main battery, O<sub>2</sub> hose, O<sub>2</sub> sensor cable, O<sub>2</sub> sensor, O<sub>2</sub> tee and operator's manual; floor stand, patient circuit support arm, 3/16" hex wrench and instructions; docking station, protective cover, power cord and operator's manual

## PTV docking stations, stands and accessories



**19270-001**

**Docking station**

Docking station, protective cover, power cord and operator manual



**13914-001**

**Basket, O<sub>2</sub> cylinder (for use w/floor stand 13912-001)**



**21380-001**

**Protective collar, docking station**



**19574-001**

**Stand, table**



**13912-001**

**Stand, floor, w/patient circuit arm**

Floor stand, patient circuit support arm, 3/16" hex wrench and instructions



**18862-001**

**Wall mount w/circuit support arm for docking station**

Wall channel, 12" support arm, mounting plate, channel cover, 16" support arm, patient circuit support arm and instructions

# ReVel ventilator accessories

## PTV docking stations, stands and accessories (continued)



**13886-001**  
Wall mount, humidifier pole



**12140-001**  
Arm, patient circuit support



**19937-001**  
Mounting plate, docking station (for use w/GCX wall mounts)



**12695-002**  
ReVel carry case

## PTV power accessories and batteries



**12103-001**  
Power adapter w/power cord, 110 VAC



**12102-001**  
Adapter, automobile power point, 12'



**13908-001**  
Battery, main lithium-ion



**10536**  
Power cord, 110 VAC, hospital grade, 10'



**13955-001**  
Charger, desktop battery w/110 VAC power cord (batteries not included)

PTV miscellaneous parts and accessories



**10699**  
Hose, O<sub>2</sub>, 15'



**13897-001**  
Cable, O<sub>2</sub> sensor



**13036-001**  
Sensor, O<sub>2</sub>



**18602-001**  
Cable, SpO<sub>2</sub>



**13037-001**  
T-Piece, reusable O<sub>2</sub> sensor



**13035-001**  
SpO<sub>2</sub> clip sensor,  
reusable adult



**13038-001**  
T-Piece, SPU O<sub>2</sub> sensor  
(10/pkg)



**19461-001**  
SpO<sub>2</sub> clip sensor,  
reusable pediatric

# ReVel ventilator accessories

## PTV miscellaneous parts and accessories (continued)



**19462-001**  
SpO<sub>2</sub> reflectance sensor



**19826-001**  
Filter, air inlet



**13063-001**  
Adapter and cable,  
VOXP data




**19833-001**  
Filter, fan (10/pkg)

Cat. No.	Size	Qty.
<b>19463-001</b>	SpO <sub>2</sub> Flexi-Form II, SPU adult	10/pkg
<b>19464-001</b>	SpO <sub>2</sub> Flexi-Form II, SPU pediatric	10/pkg
<b>21311-001</b>	SpO <sub>2</sub> Flexi-Form III, SPU adult	24/pkg

**PTV ReVel circuits**


*SPU circuits, 6'*



**14858-001**

**Patient circuit w/o PEEP w/elbow, 22mm SPU (10/pkg)**

w/1.8 m x 22 mm ID inspiratory limb, wye, elbow and 61 cm x 22 mm expiratory limb




**14859-001**

**Patient circuit w/o PEEP w/elbow, 15mm SPU (10/pkg)**

w/1.8 m x 15 mm ID inspiratory limb, wye, elbow and 61 cm x 15 mm expiratory limb


*SPU heated wire circuits for Fisher & Paykel (F&P) humidifiers, 5.5'*



**15091-102**

**Patient circuit w/o PEEP w/elbow, 22 mm SPU dual heated wire, F&P (10/pkg)**

w/1.7 m x 22 mm ID heated inspiratory limb, temp probe port, wye, elbow and 61 cm x 22 mm heated expiratory limb




**15091-104**

**Patient circuit w/o PEEP w/elbow, 15 mm SPU dual heated wire, F&P (10/pkg)**

w/1.7 m x 15 mm ID heated inspiratory limb, temp probe port, wye, elbow and 61 cm x 15 mm heated expiratory limb


*SPU heated wire circuits for ConchaTherm® or Neptune™ humidifiers, 5.5'*



**18855-001**

**Patient circuit w/o PEEP w/elbow, 22 mm SPU dual heated wire, Teleflex (10/pkg)**

w/1.7 m x 22 mm ID heated inspiratory limb, temp probe port, wye, elbow and 61 cm x 22 mm heated expiratory limb



**18856-001**

**Patient circuit w/o PEEP w/elbow, 15 mm SPU dual heated wire, teleflex (10/pkg)**


w/1.7 m x 15 mm ID heated inspiratory limb, temp probe port, wye, elbow and 61 cm x 15 mm heated expiratory limb



# ReVel ventilator accessories

## PTV ReVel circuits (continued)


Reusable patient circuits, 6'



### 14863-001

**Patient circuit w/o PEEP w/elbow, 22 mm reusable (each)**

w/1.8 m x 22 mm ID inspiratory limb, wye, elbow and 61 cm x 22 mm expiratory limb



### 14862-001

**Patient circuit w/o PEEP w/elbow, 15 mm reusable (each)**

w/1.8 m x 15 mm ID inspiratory limb, wye, elbow and 61 cm x 15 mm expiratory limb

## Circuit accessories

Cat. No.	Description	Qty.
<b>18225-001</b>	Patient wye and sense lines, reusable	
<b>18226-101</b>	Patient wye and sense lines, SPU	10/pkg

To place an order or for more information,  
please contact Customer Service at:

Tel: 1.800.754.1914

Fax: 1.763.398.8359

Email: [pulmoneticare@carefusion.com](mailto:pulmoneticare@carefusion.com)



Caution: Federal Law restricts this device to sale, distribution and use by or on the order of a physician.

CareFusion  
Yorba Linda, CA

[carefusion.com/revelcatalog](http://carefusion.com/revelcatalog)



CareFusion