

LIFEPAK[®] 15

monitor/defibrillator

Genuine accessories

LIFEPAK 15 monitor/defibrillator

Genuine accessories from **Stryker**

**Thank you for choosing Stryker as your partner
in helping save lives and improve patient care.**

You already know our reputation for quality and reliability in LIFEPAK defibrillators and AEDs. When it comes to LIFEPAK accessories, we know you won't accept anything less. Inside this catalog you'll find accessories designed to meet your needs and reliability standards. All the accessories and disposables in this catalog are designed to work with the LIFEPAK 15 monitor/defibrillator. If you have any questions about accessories or disposables, please contact your Stryker representative.



Power options



REDI-CHARGE® Base
AC power cord and adapter tray not included.
11141-000115

REDI-CHARGE LIFEPAK 15 Adapter Tray
11140-000052

AC Power Cords
11140-000015



Mobile Battery Charger
Includes AC and DC power cords, mounting bracket and operating instructions.
11577-000011



Lithium-ion Battery 5.7 Ah
5.7 amp hour, 11.1 volt, rechargeable, with fuel gauge.
21330-001176



AC Power Adapter
Includes Right Angle Cable (AC power cord not included)
11140-000098

DC Power Adapter
Includes DC cable and Right Angle Cable
11140-000074

AC Power Cord
11140-000015

Power Attachment Kit
11577-000019



Extension Cable
For AC/DC Power Adapter.
11140-000080




Replacement Right Angle Power Cable
For AC/DC Power Adapter.
11140-000081



Replacement DC Input Cable
For DC Power Adapter.
11140-000084

ECG monitoring Accessories



12-Lead ECG Cable Trunk Cable with 4-Wire Limb Leads

11111-000018 (5ft)
11111-000020 (8ft)



12-Lead ECG Cable 6-Wire Precordial Attachment

11111-000022



4-wire Cable Comb
21300-008054 (10/pack)

6-wire Cable Comb
21330-008055 (10/pack)



3-Wire ECG Cable

11110-000029



5-Wire ECG Cable

Right-angle connector, 4-wire limb plus one chest lead, labeled "V1" on the LIFEPAK 15 monitor reports.

11110-000066



Box of Strip Chart Recorder Paper

100mm x 22m
11240-00003 2 (2 rolls per box)



LIFE•PATCH® ECG Electrodes

Adult, pregelled
11100-000001 (3/pack)



LIFE•PATCH® ECG Electrodes

Adult, pregelled
11100-000002 (4/pack)

Therapy delivery accessories

Hard paddles and electrode gel



Standard Hard Paddles
11130-000061 (1 pair)



Pediatric Paddle, External
(Two required); slips onto standard adult hard paddle.
11133-000007 (1 adapter)



SIGNAGEL® Electrode Gel
For use with hard paddles. Highly conductive, multipurpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback and EMG.
21300-005847 (8.5 oz)



Internal Paddles (Requires Internal Paddle Handles and Internal Paddles Adapter Cable)

1 in. size 11131-000010 (1 pair, 6.25 in. shaft)	2.5 in. size 11131-000013 (1 pair, 5.5 in. shaft) 11131-000019 (1 pair, 8.5 in. shaft)
1.5 in. size 11131-000011 (1 pair, 6 in. shaft) 11131-000021 (1 pair, 9 in. shaft) 11131-000024 (1 pair, 14 in. shaft)	3.5 in. size 11131-000014 (1 pair, 5 in. shaft) 11131-000023 (1 pair, 8 in. shaft)
2 in. size 11131-000012 (1 pair, 5.75 in. shaft) 11131-000022 (1 pair, 8.75 in. shaft)	



Internal Paddle Handles with Discharge Control
(For use with the Internal Paddles Adapter Cable)
11131-000001 (1 pair)



Internal Paddles Adapter Cable
(For use with Internal Paddle Handles)
11998-000326

Therapy delivery accessories

EDGE System™ Electrodes for pacing/defibrillation/ECG with QUIK-COMBO® Connector

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector

24 in. leadwire length
11996-000091



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector

24 in. leadwire length
11996-000090



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK® Preconnect System

42 in. leadwire length
11996-000017



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only with manual monitor/defibrillators; 12 month minimum shelf life at time of shipment 24 in. leadwire length.
11996-000093



QUIK-COMBO Therapy Cable

With convenient TRUE-LOCK® Cable Connector. Length is approximately 8 ft. For use with LIFEPAK 15 monitor/defibrillator.
11113-000004 (8ft)

NIBP monitoring accessories

NIBP hoses



NIBP Tubing, Straight

21300-008159 (6ft)
21300-008147 (9ft)
21300-008146 (12ft)



NIBP Tubing, Coiled

21300-008148 (2-9ft)

NIBP cuffs



Reusable Cuff

X-Large Adult 35 - 44 cm 11160-000019	Pediatric 13 - 20 cm 11160-000013
Large Adult 32 - 42 cm 11160-000017	Infant 8 - 14 cm 11160-000011
Adult 2 26 - 35 cm 11160-000015	



Single Patient Use Cuff

X-Large Adult 35 - 44 cm 11160-000020	Pediatric 13 - 20 cm 11160-000014
Large Adult 32 - 42 cm 11160-000018	Infant 8 - 14 cm 11160-000012
Adult 26 - 35 cm 11160-000016	

Pulse oximetry monitoring accessories

Masimo SET® RC patient cables



RC Patient Cable
For use with M-LNCS and
Rainbow Patient Sensors
11171-000037 (4ft)
11171-000038 (12ft)



RC EMS Patient Cable
11171-000082 (4 ft)

Masimo SET RC patient cable compatible SpO₂ sensors



M-LNCS Reusable Sensor
11171-000046 (Ad)
11171-000047 (Ped)




**M-LNCS Adhesive Sensors
(20/box)**
11171-000039 (Ad)
11171-000040 (Ped)



**M-LNCS Adhesive Sensors
(20/box)**
11171-000043 (Neo/Pt)
11171-000042 (Neo/Ad)
11171-000041 (Inf)

Masimo SET RC patient cable compatible Rainbow® SpO₂, SpCO, SpMet sensors



Rainbow Reusable Sensor
11171-000049 (Ad)
11171-000050 (Ped)



**Rainbow Adhesive Sensor
(10/box)**
11996-000342 (Inf)
11996-000341 (Neo/Ad)



**Rainbow Adhesive Sensor
(10/box)**
11996-000339 (Ad)
11996-000340 (Ped)

Pulse oximetry monitoring accessories

Masimo SET LNC patient cables



Red LNC Patient Cable

For Use with LNCS Patient Sensors

11996-000323 (4ft)

11996-000324 (10ft)

11996-000325 (14ft)

Masimo SET LNC patient cable compatible SpO₂ sensors



LNCS® Reusable Sensor

11171-000017 (Ad)

11171-000018 (Ped)



LNCS Reusable Soft Sensor

11171-000052 (Ad)



LNCS Adhesive Sensor (20/box)

11171-000019 (Ad)

11171-000020 (Ped)



LNCS Adhesive Sensor (20/box)

11171-000029 (Neo/Pt)

11171-000028 (Neo/Ad)

11171-000031 (Inf)

Pulse oximetry monitoring accessories

Direct connect SpO₂ only patient sensors



Adult Reusable Direct Connect Sensor

11996-000331 (3ft)
11996-000332 (12ft)



Pediatric Reusable Direct Connect Sensor

11996-000333 (3ft)
11996-000334 (12ft)



Adult Reusable Soft Direct Connect Sensor

11171-000053 (8ft)

Direct connect Rainbow SpO₂, SpCO, SpMet patient sensors



Adult Rainbow Direct Connect Reusable Sensor


11996-000335 (3ft)
11171-000032 (8ft)
11996-000336 (12ft)



Pediatric Rainbow Direct Connect Reusable Sensor

11996-000337 (3ft)
11171-000033 (8ft)
11996-000338 (12ft)

Additional Masimo accessories



Reusable Ambient Light Shield

11171-000054 (5/bag)



Disposable Ambient Light Shield

11171-000055 (10/bag)

Pulse oximetry monitoring accessories

Masimo to Nellcor™ adapter



Red MNC Cable

Connects LIFEPAK 15 to
Nellcor patient sensor
11996-000365 (4ft)
11996-000366 (10ft)

Masimo to Nellcor adapter compatible SpO₂ sensors



Nellcor Reusable SpO₂ Sensor, Adult, DS100A

11996-000060 (Adult - Ref# DS100A)



Nellcor Reusable Multisite SpO₂ Sensor, D-YS

11996-000106 (≥ 1 kg - Ref# D-YS)



Nellcor Reusable two-piece SpO₂ sensors, OXI-A/N, OXI-P/I

Includes 50 disposable adhesive sensors.
11996-000061 (Ad/Neo)
11996-000062 (Ped/Inf)

Nellcor Single-Patient Use Adhesive Sensor Wraps (100/pk)

Not Available in Canada.
11996-000048 (Ad/Neo - Ref# ADH A/N)
11996-000049 (Ped/Inf - Ref# ADH P/I)

Nellcor Single-Patient Use Foam Wrap

Not Available in Canada.
11996-000110 (Ad/Neo - Ref# FOAM A/N)
11996-000108 (Ped/Inf - Ref# FOAM P/I)

Temperature monitoring

Temperature cables



Temperature Adapter Cable

11140-000079 (10ft)

11140-000078 (5ft)

Temperature sensors



Esophageal-Rectal Temperature Sensor

11996-000360 (9FR, 20/box)



Skin Temperature Sensor

11996-000359 (20/box)

End-Tidal CO₂ (EtCO₂) monitoring accessories

Microstream® filterlines for intubated patients

Single patient use



FilterLine® SET

Key Applications: OR, EMS, ED, Rapid Response Teams, Transport.

Adult/Pediatric

11996-000081 (25/pk, 2 m)

11996-000164 (25/pk, 4 m)



FilterLine H SET

Key Applications: Critical care, Humidified Environments.

Adult/Pediatric

11996-000080 (25/pk, 2 m)



FilterLine H SET

Key Applications: Critical care, Humidified Environments.

Infant/Neonatal

11996-000001 (25/pk, 2 m)

Microstream non-intubated filterlines

Single patient use



Smart CapnoLine® Plus

Key Applications: Procedural Sedation, Upper GI Procedures, MAC, RMS, BD, Rapid Response Teams.

Adult with O₂

11996-000163 (25/pk, 2 m)

11996-000167 (100/pk, 2m)

11996-000165 (25/pk, 4 m)

Adult without O₂

11996-000162 (25/pk, 2 m)

11996-000166 (100/pk, 2 m)



Smart CapnoLine

Key Applications: Procedural Sedation, SMS, SD, Rapid Response Teams.

Pediatric with O₂

11996-000128 (25/pk, 2 m)

Pediatric without O₂

11996-000120 (25/pk, 2 m)

LIFEPAK 15 monitor/defibrillator

Cases & mounting options



Standard Carrying Case

Includes right pouch and left pouch, for use with LIFEPAK 15 monitor/defibrillator.
11577-000002

Top Pouch

Storage for sensors and electrodes; insert in place of standard paddles.
11220-000028

Shoulder Strap

For use with LIFEPAK 15 monitor/defibrillator.
11577-000001

Back Pouch

Ideal for additional accessory storage. For use with LIFEPAK 15 monitor/defibrillator.
11260-000039



Bed Connector

For use in hospital only.
11996-000374

Communication accessories



LIFEPAK Monitor to PC Cable

For connecting LIFEPAK 12 or LIFEPAK 15 monitor/defibrillator to PC.
11230-000020 (Serial)
11996-000369 (USB)



Titan III Wireless Gateway

For transmitting data from LIFEPAK 15 to the LIFENET® System or CODE-STAT™ data review software. Requires existing wireless network.
21996-000109 (Wireless Gateway)
11996-000479 (Wireless/Cellular/Audio Gateway) - For use with Stryker Verizon Data Plan
11996-000484 (Wireless/Cellular/Audio Gateway) - For use with Stryker AT&T Data Plan
11996-000480 (Wireless/Cellular/Audio Gateway) - For use with customer data plan



4G Gateway

For transmitting data from LIFEPAK 15 to the LIFENET System or CODE-STAT data review software. Requires data plan.
11996-000476 (AT&T Gateway)
11996-000475 (AT&T Gateway - for use with Stryker Data Plan)
11996-000474 (Verizon Gateway)
11996-000471 (Verizon Gateway - for use with Stryker Data Plan)

Stryker Data Plans

99428-000305 Verizon Data Plan 1-Year
99428-000304 AT&T Data Plan 1-Year

Testers & training materials




12-Lead Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias, and bradycardias, as well as ST segment and T wave abnormalities. For use with LIFEPAK devices with a 12-lead ECG feature.
11996-000311 (12-Lead)

3-Lead Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe, interactive training. Select from 17 ECG rhythms including: fibrillation, tachycardias and bradycardias.
11996-000310



Defibrillator Checker

Tests integrity of energy delivery through standard hard paddles. Neon light indicates energy has been delivered.
11998-000060



Test Load

Use to perform therapy cable performance checks. Connects to the QUIK-COMBO therapy cable on the defibrillator.
21330-001365 (English only)
21330-001367 (English, French for Canada)

Training tools



Operating Instructions: LIFEPAK 15 Monitor/Defibrillator

Also available as free online download at <https://www.strykeremergencycare.com/disclosure-and-safety-information/>
26500-002408

**For further information, please contact Stryker at 800 442 1142 (U.S.),
800 668 8323 (Canada) or visit our website at strykeremergencycare.com**

Emergency Care

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CODE-STAT, EDGE System, LIFEPAK, REDI-PAK, QUIK-COMBO, REDI-CHARGE, Stryker, TRUE-LOCK, Masimo, the Radical logo, Rainbow and SET are registered trademarks of Masimo Corporation. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

GDR 3300304_N
Copyright © 2020 Stryker

Distributed by:

Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
strykeremergencycare.co

Distributed by:

Stryker Canada
2 Medicorum Place
Waterdown, Ontario
L8B 1W2
Canada
Toll free 800 668 8323