

LIFEPAK[®] 20e

defibrillator/monitor with
CodeManagement Module[®]

Genuine accessories

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module





Genuine accessories from **Stryker**

Ensure the safety of your staff and patients.

Only Stryker accessories and disposables have been thoroughly tested with LIFEPAK products.

Stick with the name you trust and the company that stands behind it. To order, contact your Stryker representative.

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

Power options



AC Power Cord

Powers LIFEPAK 20e defibrillator/
monitor and recharges the
internal battery.

11140-000015



Rechargeable Internal Battery

For use as AC power backup.

11141-000112 (LIFEPAK 20e Lithium-ion)



Lithium-ion Battery for CodeManagement Module

Rechargeable battery for the CodeManagement Module; for
use as AC power backup.

11141-000162 (Lithium-ion)

ECG monitoring accessories



5-Wire ECG Cable
11110-000066



3-Wire ECG Cable
11110-000029



LIFE-PATCH® ECG Electrodes
Adult, pregelled.
11100-000001 (3/pack)
11100-000002 (4/pack)

Cleartrace® ECG Electrodes
Adult, clear tape.
Available in the EU.
1700-003 (3/pack)



50mm x 22mm ECG Printer Paper
11240-000031 (3 rolls/box)



CodeManagement Module for use with LIFEPAK 20e Defibrillator/Monitor - Wireless
Connects to the LIFEPAK 20e defibrillator to add wireless connectivity.
11150-000018

CodeManagement Module for use with LIFEPAK 20e Defibrillator/Monitor - Wireless & Capnography
Connects to the LIFEPAK 20e defibrillator to add wireless connectivity and capnography monitoring.
11150-000019

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

Therapy delivery accessories

Hard paddles and electrode gel



Standard Adult Detachable Hard Paddles

Adult paddle plates are removable to reveal pediatric paddles.
11130-000037



SIGNAGEL® Electrode Gel

For use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel.
21300-005847 (8.5 oz)



Internal Paddle Handles with Discharge Control

For use with LIFEPAK 15, 20e defibrillator/monitors.
11131-000001 (1 pair)



Internal Paddles (Requires Internal Paddle Handles and Internal Paddles Adapter Cable)

1 in. size
11131-000010 (1 pair, 6.25 in. shaft)

1.5 in. size
11131-000011 (1 pair, 6 in. shaft)
11131-000021 (1 pair, 9 in. shaft)
11131-000024 (1 pair, 14 in. shaft)

2 in. size
11131-000012 (1 pair, 5.75 in. shaft)
11131-000022 (1 pair, 8.75 in. shaft)

2.5 in. size
11131-000013 (1 pair, 5.5 in. shaft)
11131-000019 (1 pair, 8.5 in. shaft)

3.5 in. size
11131-000014 (1 pair, 5 in. shaft)
11131-000023 (1 pair, 8 in. shaft)

Therapy delivery accessories

EDGE System™ Electrodes for pacing/defibrillation/ECG with QUIK-COMBO® Connector

18 month minimum shelf life unless specified



EDGE System Electrodes with QUIK-COMBO Connector
11996-000091 (24" leadwire length)



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector
11996-000090 (24" leadwire length)



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only with manual monitor/defibrillators; 12 month minimum shelf life at time of shipment.

11996-000093 (24" leadwire length)



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK® Preconnect System
11996-000017 (42" leadwire length)



QUIK-COMBO Therapy Cable
11110-000004

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

End-Tidal CO₂ (EtCO₂) monitoring accessories

Filterlines for intubated patients

Single patient use



FilterLine® SET

Key Applications: OR, EMS, ED, Rapid Response Teams, Transport.

Adult/Pediatric

11996-000081 (25/pk, 200 cm)

11996-000164 (25/pk, 400 cm)



FilterLine H SET

Key Applications: Critical care, Humidified Environments.

Adult/Pediatric

11996-000080 (25/pk, 200 cm)

Infant/Neonatal

11996-000001 (25/pk, 200 cm)

Non-intubated filterlines

Single patient use



Smart CapnoLine® Plus

Key applications: Procedural sedation, Lower GI procedures, MAC, EMS, ED, Rapid Response Teams.

Adult with O₂

11996-000163 (25/pk, 200 cm)

11996-000167 (100/pk, 200 cm)

11996-000165 (25/pk, 400 cm)

Adult without O₂

11996-000162 (25/pk, 200 cm)

11996-000166 (100/pk, 200 cm)



Smart CapnoLine

Key applications: Procedural sedation, EMS, ED, Rapid Response Teams.

Pediatric with O₂

11996-000128 (25/pk, 200 cm)

Pediatric without O₂

11996-000120 (25/pk, 200 cm)

SpO₂ monitoring accessories

Masimo LNC® reusable patient cables



LNCS Patient Cable

11171-000024 (4 ft)
 11171-000016 (10 ft)
 11171-000025 (14 ft)
 11171-000027 (4 ft extension)

Masimo LNC compatible SpO₂ only sensors



LNCS® Adhesive Sensor

Single Patient Use. (20/box)
 11171-000019 (Ad)
 11171-000020 (Ped)
 11171-000029 (Neo/Pt)
 11171-000028 (Neo/Ad)
 11171-000031 (Inf)



LNCS® Reusable Sensor

11171-000017 (Ad)
 11171-000018 (Ped)



LNCS Reusable Soft Sensor

11171-000052 (Ad)

Masimo reusable direct connect sensors and cables



Direct Connect Reusable Soft Sensor

11171-000051 (Ad)

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

SpO₂ monitoring accessories

Additional accessories and adapter cables



**Reusable Ambient
Light Shield**
11171-000054 (5/pack)



**Disposable Ambient
Light Shield**
11171-000055 (10/pack)

SpO₂ monitoring accessories

Masimo to Nellcor™ adapter



MNC-1 Adapter Cable

Allows specially configured LIFEPAK 20e defibrillator/monitor with Masimo SpO₂ to connect to Nellcor sensors.

11996-000183 (10 ft)

11996-000198 (4 ft)

Masimo to Nellcor adapter compatible SpO₂ sensors



DURASENSOR Reusable Clip

11996-000060 (Ad)

DURA-Y Multisite Reusable Sensor

11996-000106 (> 1 kg)



Oxiband Reusable Sensor

Includes 50 disposable adhesive sensors.

11996-000061 (Ad/Neo)

11996-000062 (Ped/Inf)



Disposable Adhesive Bandage Wrap (100/pk)

Not Available in Canada.

11996-000048 (Ad/Neo)

11996-000049 (Ped/Inf)



Oxisensor II Adhesive Sensors

(24/box)

11996-000113 (Ad, 18 in.)

11996-000114 (Ad, 36 in.)

11996-000116 (Ped, 18 in.)



Oxisensor II Adhesive Sensors

(24/box)

11996-000115 (Inf, 18 in.)

11996-000117 (Neo/Ad 18 in.)

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

Carrying cases and mounting options



Basic Carrying Case/Accessories Organizer

Includes set of right and left pouches for additional accessories.

11260-000018



Carrying Case for LIFEPAK 20e Defibrillator with Module

Includes a set of right and left pouches for additional storage of accessories. For use with a LIFEPAK 20e device with CodeManagement Module.

11260-000045



Top Pouch

For additional storage: insert in place of standard paddles.

11260-000043

Shoulder Strap for Basic Carrying Case

Includes mounting brackets and hardware.

11260-000041



Docking Station

Attaches to crash cart to lock LIFEPAK 20e defibrillator/monitor safely in place.

21330-000996

Right Side Accessory Pouch

Attaches to LIFEPAK 20e defibrillator/monitor and holds items such as QUIK-COMBO therapy cable, electrodes and SIGNAGEL electrode gel.

11260-000016

Training tools and testers




Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias and bradycardias, as well as ST segment and T wave abnormalities.

11996-000311 (12-Lead)

11996-000310 (3-Lead)



Defibrillator Checker

Tests integrity of energy delivery through the standard hard paddles. Neon light indicates energy has been delivered.

11998-000060



QUIK-COMBO Test Plug

Test plug connects to QUIK-COMBO therapy cable to test cable integrity.

11113-000002

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

Literature, communication and third party accessories



Inservice DVD

26500-003545 (English - DVD)

Service Manual

26500-003543 (English - CD - LIFEPAK 20e
with CodeManagement Module)

Operating Instructions:

LIFEPAK 20e Defibrillator/Monitor
with CodeManagement Module
26500-002570

Communication accessories



Serial Port Cable

For connecting LIFEPAK 20e
defibrillator/monitor to PC.
11230-000018 (6 ft)

Transport Configuration Cable

For connecting LIFEPAK 20e to
LIFEPAK 20e defibrillator/monitor.
11230-000019 (6 ft)



**For further information, please contact Stryker at 800 442 1142 (U.S.),
800 668 8323 (Canada) or visit our website at strykeremergencycare.com**

Emergency Care

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CodeManagement, EDGE System, LIFEPAK, REDI-PAK, QUIK-COMBO, Stryker, Masimo, the Radical logo, Rainbow and SET are registered trademarks of Masimo Corporation. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

GDR 3300292_J

Copyright © 2020 Stryker

Distributed by:

Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
strykeremergencycare.com

Distributed in Canada by:

Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
strykeremergencycare.com