



## SPECIFICATIONS

Product Description: **PLASTIC DEFIBRILLATOR WALL CASE W/ ALARM**  
 Part Number: **4400-DA**  
 Style: **SMALL, SINGLE DOOR**

### GENERAL DESCRIPTION:

White case made of corrosion resistant ABS plastic and high visibility defibrillator graphics. Special features include 6" x 10" (15.2 cm x 25.4 cm) window and right-hand door opening. Three optional slide-in shelves. Comes with 9-volt battery powered alarm that goes off when the door is opened.

### CONSTRUCTION:

Door: Single hinged, riveted, white color  
 Seal: Buna-N Rubber  
 Material: Acrylonitrile Butadiene Styrene (ABS) with ultraviolet inhibitors  
 Case: Black color

### WINDOW:

Material: Cast Acrylic, Clear  
 Size: 8" x 12" (20.3 cm x 30.4 cm)

### ALARM:

Power supply: 9VDC (Battery)  
 Loudness: 130 Decibel

### SHELF: (Sold Separately)

Slots: Case has three slots for shelf position  
 Material: Clear cast acrylic  
 Part No: 4511-00 sold by each  
 Size: 12<sup>3</sup>/<sub>4</sub>" x 12<sup>1</sup>/<sub>4</sub>" x 6<sup>1</sup>/<sub>4</sub>" (32.3 cm x 31.1 cm x 15.8 cm)



### DIMENSIONS:

Exterior			Interior			Weight	
Width	Depth	Height	Width	Depth	Height	Shipping	Product
14" (35.5 cm)	9 <sup>1</sup> / <sub>2</sub> " (23.4 cm)	18" (45.7 cm)	12" (30.4 cm)	7 <sup>1</sup> / <sub>2</sub> " (19.05 cm)	16" (40.6 cm)	6 lbs. (2.7 kg)	4 lbs. (1.8 kg)