



## SPECIFICATIONS

Product Description: **PLASTIC DEFIBRILLATOR WALL CASE W/STROBE**  
 Part Number: **4400-DS**  
 Style: **SMALL, SINGLE DOOR**

### GENERAL DESCRIPTION:

White case made of corrosion resistant ABS plastic and high visibility defibrillator graphics. Special features include 6" (15.2 cm) x 10" (25.4 cm) window and right hand door opening. Three optional slide in shelves. Comes with a strobe and alarm that goes off when door is opened.

### CONSTRUCTION:

Door: Single hinged, riveted, white color  
 Seal: Buna-N Rubber  
 Material: Acrylonitrile Butadiene Styrene (ABS) with ultraviolet inhibitors  
 Case: Black color

### WINDOW:

Material: Cast Acrylic, Clear  
 Size: 8" (20.3 cm) x 12" (30.4 cm)

### STROBE:

Power supply: 9 V DC Battery  
 Features: Fully sealed, water-proof

### ALARM:

Power supply: 9 V DC Battery  
 Loudness: 120 Decibel

### SHELF: (Sold Separately)

Slots: Case has three slots for shelf position  
 Material: Clear cast acrylic  
 Part No: 4511-00 sold by each  
 Size: 12¾" (32.3 cm) X 12¼" (31.1 cm) X 6¼" (15.8 cm)



### DIMENSIONS:

Exterior			Interior			Weight	
Width In (cm)	Depth In (cm)	Height In (cm)	Width In (cm)	Depth In (cm)	Height In (cm)	Shipping Lbs (Kg)	Product Lbs (Kg)
14" (35.5 cm)	9½" (23.4 cm)	18" (45.7 cm)	12" (30.4 cm)	7½" (19.05 cm)	16" (40.6 cm)	6 lbs (2.7 kg)	4 lbs (1.8 kg)