



V-VAC Starter Kit for Patient Use (shown above) includes:
 Handle, 2 Replacement Cartridges, Short Suction Catheter w/Adapter Tip, Double Male Connector, and Directions for Use.

Ordering information:

- 98 50 00 V-VAC Starter Kit
- 98 51 00 V-VAC Training Kit
- 98 50 01 Replacement Cartridge
- 98 52 00 Carrying Bag
- 98 50 30 Replacement Handle

V-VAC Training Kit for Practical Use includes:

Handle, 1 Training Cartridge with plug for emptying and cleaning, Catheter w/Adapter Tip and Double Male Connector.*

*Not illustrated.



A The V-VAC unit fits easily into many emergency kits. A separate carry bag is also available.

B The soft tip of the intake valve has deep cutouts to prevent tube occlusion and soft tissue injury. The tri-cut valve resists back flow of suctioned matter.

For additional information or to order, contact your local Laerdal Company or Distributor



Norway: LAERDAL MEDICAL AS
 PO Box 377, N-4002 Stavanger
 Tel: +47 51 51 17 00, Fax +47 51 52 35 57
 E-mail: laerdal.norway@laerdal.no

Australia: LAERDAL Pty Ltd. (ABN 47 003 817 490)
 23 Edward Street, Oakleigh, Victoria 3166
 Tel: +61 3 9569 4055, Toll free 1800 331 565
 Fax +61 3 9569 4028, Toll free 1800 635 835
 E-mail: customerservice@laerdal.com.au

Benelux: "The Laerdal BeNeLux House"
 Dijkstraat 54, NL-5554 PS VALKENSWAARD
 Tel: +31 40 208 58 00, Fax +31 40 208 58 02
 E-mail: info@laerdal.nl

Benelux: Belgium and Luxembourg
 LAERDAL BeNeLux nv
 Mechelsesteenweg 277, B-1800 VILVOORDE
 Tel: +32 2 253 36 96, Fax + 32 2 253 36 90
 E-mail: info@laerdal.be

Canada and Latin America:
 LAERDAL MEDICAL CANADA LTD.
 151 Nashdene Rd., Unit #45
 Toronto, ON, Canada, M1V 4C3
 Tel: +1 (416) 298-9600,
 Toll free 888/LAERDAL (523-7325)
 ou en français (800) 567-9987
 Fax +1 (416) 298-8016
 E-mail: savelves@laerdal.ca

Denmark: LAERDAL DANMARK
 Esplanaden 18, 1
 DK-1263 Copenhagen K
 Tel: +45 80 333 112, Fax: +45 80 333 555
 E-mail: laerdal.denmark@laerdal.no

Eastern Europe: LAERDAL C.E. EUROPE
 A-1060 Vienna, Mariahilferstr. 1/d/1/3
 Tel-Fax +43 15877140, Tel: +43 15815927
 E-mail: mstam@aon.at

Far East: LAERDAL Singapore Pte Ltd
 No 1 Marine Parade Central
 #13-05, Parkway Centre, Singapore 449408
 Tel: +65 63464259, Fax: +65 63467523
 E-mail: admin@laerdal.com.sg

Finland: LAERDAL OY
 Vattunienranta 2, 00210 HELSINKI
 Puhelin +358 (0)9-612 99 80,
 Fax +358 (0)9-692 77 99
 E-mail: laerdal.finland@laerdal.no

France: LAERDAL MEDICAL FRANCE
 1 rue des Vergers - Bâtiment n° 5 - 69578 Limonest Cedex
 Tel: +33 (0)4 72 52 02 52, Fax: +33 (0)4 78 35 38 45
 E-mail: laerdal.france@laerdal.no

Germany: LAERDAL MEDICAL DEUTSCHLAND GmbH
 Am Lohrlefeld 56, 81249 München
 Tel: +49 (0)89 / 864 954-0, Fax +49 (0)89 / 864 34 84
 E-mail: info@laerdal.de

Italy: LAERDAL ITALIA s.r.l.
 Via Piero Gobetti, 5/22, int.Z, 40129 Bologna
 Tel: +39 051-355587, Fax: +39 051-355598
 E-mail: laerdal.italy@laerdal.no

Japan: LAERDAL MEDICAL JAPAN K.K.
 Shinjuku IS Bldg., 8F
 2-13-12, Shinjuku, Shinjuku-ku, Tokyo 160-0022
 Tel: 81-3-5919-2121, Fax: 81-3-5919-2122
 E-mail: laerdal.japan@laerdal.co.jp

Malaysia: LAERDAL HOSPILINE Sdn Bhd
 12 Jalan Titar 33/25, Section 33, Shah Alam Technology Park,
 40400 Shah Alam, Selangor Darul Ehsan,
 Tel: +60 (3) 5122-7002, Fax +60 (3) 5122-5650
 E-mail: customerservice@laerdal.com.my

New Zealand: LAERDAL NEW ZEALAND Ltd
 (GST Number 79-403-636)
 Po Box 36339, Merivale, Christchurch
 Tel: +61 3 9564 4066, Toll free 0800 523 732
 Fax +61 3 9563 3368, Toll free 0800 528 852
 E-mail: customerservice@laerdal.co.nz

Portugal: LAERDAL ESPANA, Sucursal em Portugal
 Rua da Bela Vista à Graca, no.31 A, Escritorio 9,
 1170-054 Lisboa
 Tel: (+351)-21-816-6567, Fax (+351)-21-816-6568,
 Mob. (+351)-918765367

Spain: LAERDAL ESPANA, S.L.
 Manuel Tovar 19, 28034 Madrid
 Tel: (+34) 902-291110, Fax: (+34) (91)-7291125
 E-mail: laerdal.spain@laerdal.no

Sweden: LAERDAL MEDICAL AB
 Box 2064, Paviljongvägen 5, 132 02 Saltsjö-Boo
 Tel: +46 (0)8-55614610, Fax +46 (0)8-55614619
 E-mail: kontakt@laerdal.se

United Kingdom: LAERDAL MEDICAL LTD.
 Laerdal House, Goodmead Road, Orpington,
 Kent BR6 0HX
 Tel: +44 (0)1689 876634, Fax: +44 (0)1689 873800
 E-mail: customerservice@laerdal.co.uk

USA:
 LAERDAL MEDICAL CORPORATION
 167 Myers Corners Road, PO Box 1840
 Wappingers Falls, New York 12590-8840
 Tel: (800) 431-1055, +1 (845) 297-7770
 Fax: (800) 227-1143, +1 (845) 298-4545
 E-mail: customerservice@laerdal.com

LAERDAL TEXAS
 PO Box 38, 226 FM 116 South Gatesville,
 Texas 76528-0038
 Tel: (800) 433-5539, (254) 865-7221
 Fax: (254) 865-8011
 E-mail: customerservice@laerdal.com

Get a good grip on handling even the toughest suction scenario



V-VAC Manual Suction Unit

Ideal for First Responder emergency kits and as back-up suction for hospitals, clinics and emergency responders

- All-in-one system adjusts for high and low vacuum suctioning
- Efficient: exceeds 380 mmHg suction with a peak air flow of over 70 litres per minute
- Unique, soft, no clog suction tip
- Catheter connection allows for versatile use
- Disposable, one-piece cartridge eliminates the need for cleaning contaminated parts
- Reusable handle is easy to clean for repeated use

