

Evacuation Chair

for Public Access

The complete evacuation solution

Stryker's Evacuation Chair is essential for the evacuation of disabled persons weighing up to 500 lb (227 kg), from multilevel facilities in emergency situations. The patented Stair-TREAD™ tracks allow easy evacuation down stairs by a single operator, for passengers weighing up to 200 lb (90 kg). The extendable handle allows the operator to have a forward-facing, upright posture for balance, line-of-sight and ergonomic safety. Foot end lift handles also extend, allowing two people to lift an evacuee over obstacles or up stairs.

Features and benefits

- Easy set-up and instruction pictograms on seatback
- Innovative Stair-TREAD system allows easily controlled stair descent
- Color-coded controls for faster operation
- Extendable head end and foot end lift handles for ergonomic lifting
- Extendable upper control handle for upright operator posture and unobstructed line of sight
- Molded hand grips fit hands for better ergonomics
- Positive action locking mechanism for passenger security
- Dual wheel locks for more secure transfers
- Oversized rear wheels with sealed bearings for smoother transports
- Compact storage size conserves space in corridors and closets
- Lightweight, durable, aluminum construction for low effort transports
- Power washable for easy cleaning
- Grease-free maintenance for ease of care
- Restraints with Fastex® buckles for passenger security
- Head support allows evacuation of passengers with compromised head control
- Two-piece molded ABS seat for greater passenger comfort



Optional accessories

Product description	Part number
Aggressive Stair-TREAD belt*	6253-070-000
Foot rest	6254-029-000
Wall bracket	6254-050-000
Storage cabinet	6254-002-000
Chair cover	6254-055-000
Storage cabinet alarm	6254-045-000
O₂ bottle holder	6250-140-000

*The Aggressive Stair-TREAD belt is designed to improve the grip of the Stair-TREAD system on steps that have a smooth surface with a rounded edge.

Evacuation Chair

Key features

- 1 Restraints**
with Fastex buckles for passenger security
- 2 Upper control handle**
extends to be ergonomically positioned for optimal operator control
- 3 Head end lift handles**
extend for upright lifting posture
- 4 Patented Stair-TREAD system**
controls descent in proportion to the passenger's weight
- 5 Wheels**
provide maneuverability in narrow passageways
- 6 Foldaway foot support**
provides additional passenger security and stability
- 7 Foot end lift handles**
extend for ergonomic lifting
- 8 Contoured seat**
for comfort



Warranty

15 years

Specifications

Evacuation Chair

Part number	6254-000-000
Height	45 in (114 cm)
Width	20.5 in (52 cm)
Depth¹	28 in (71 cm)
Folded depth	11 in (28 cm)
Weight	34 lb (15.5 kg)
Max. weight capacity^{2,3}	500 lb (227 kg)

Foot rest option

Part number	6254-029-000
Step height from floor	4.5 in (11.4 cm)
Floor clearance	2.25 in (5.7 cm)
Width	12.25 in (31.1 cm)
Depth	3.5 in (8.9 cm)

Storage cabinet

Part number	6254-002-000
Height	47 in (119.4 cm)
Width	24 in (61 cm)
Depth	10.25 in (26 cm)
Weight	42 lb (19 kg)

1. Depth dimensions are measured with extendable handles retracted.

2. Operators must consider the weight of the passenger, equipment and accessories when determining the total weight on the product.

3. Maximum weight capacity is total weight distributed in accordance with basic human anatomy.

For further information, please contact local Stryker sales representative or visit our website at stryker.com

The information presented is intended to demonstrate Stryker's product offerings. Refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Stair-TREAD, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

Mkt Lit-2081 07 OCT 2020 REV A
Copyright © 2020 Stryker

Manufactured by:

Stryker
3800 E. Centre Avenue
Portage, MI 49002 U.S.A.
Toll free 800 784 4336
stryker.com

Distributed in Canada by:

Stryker Canada
2 Medicorum Place
Waterdown, Ontario
L8B 1W2
Canada
Toll free 800 668 8323